

BeneHeart R12

Electrocardiograph



Technical Specifications

Physical

Height:	128 mm
Width:	365 mm
Depth:	305 mm
Weight:	4.8 Kg

Processing

ECG sampling rate:	1000 samples/second/channel
Pacer detection sampling rate:	16,000 samples/second/channel
ECG amplifier:	DC-coupled
Acquisition mode:	Pre-acquisition or post-acquisition, provide 10 seconds of instantaneous ECG acquisition
Dynamic range:	AC Differential ± 10 mV, DC offset ± 600 mV
Resolution:	1 μ V/LSB
Frequency response:	-3 dB @ 0.05 to 150 Hz
Baseline drift filter:	0.05 Hz, Baseline Drift Removal (BDR)
High cut-off frequency:	20 Hz, 35 Hz, 150 Hz
AC filter:	50/60 Hz
Common mode rejection:	≥ 110 dB (with AC filter switched off)
ADC:	24 bit
Input impedance:	≥ 50 M Ω @ 10 Hz, defibrillator protected
Time Constant:	≥ 3.2 s
Noise Level:	≤ 15 μ V
Patient leakage:	< 10 μ A
Heart rate meter:	30 to 300 BPM $\pm 10\%$ or ± 5 BPM, whichever is greater.
Sensitivity/gain:	2.5, 5, 10, 20, L=10 C=5, L=20 C=10 mm/mV, Auto
Boot time:	≤ 7 s

Display

Color:	24 bit Color display
Display type:	8 inch diagonal, TFT LCD with LED graphics backlit
Display resolution:	800x480
Display data:	Patient ID, Patient name, gender, age, heart rate, clock, battery power indicator, waveforms, lead labels, speed, gain, filter settings, warning messages, information messages, network status, USB

Power

Power supply:	AC input (without external power adaptor) or battery operation
---------------	--

AC Power

Input voltage:	100 to 240 VAC $\pm 10\%$
Input power:	100 VA
AC frequency:	50/60 Hz ± 3 Hz

Battery

Battery type:	Rechargeable Lithium ion battery, 4500 mAh
Charge time:	≤ 6 h to 90% and ≤ 7 h to 100% with device off
Battery Capacity:	More than 3.5 hours of continuous operation or 400 ECGs recording

Writer	
Writer technology:	Thermal dot array
Writer speed:	5, 12.5, 25 & 50 mm/s
Number of traces:	12 leads
Writer speed accuracy:	±5%
Writer amplitude accuracy:	±5%
Writer resolution:	Horizontal 40 dots/mm @ 25 mm/s, Vertical 8 dots/mm
Paper type:	Thermal Z-fold A4 paper (210 mm x 295 mm)
	US Letter 8.5x11 in (215 mm x 280 mm)
Software	
Measurement and interpretation:	Supports measurement and interpretation with the University of Glasgow ECG analysis algorithm
Resting ECG mode:	Records and prints 12-lead resting ECG with 10-second duration as a standard feature
Supported patient information:	Patient ID, secondary ID, name, age, date of birth, gender, race, medication, class, V3 electrode Placement.
Internal storage:	800 ECGs in internal memory
ECG Storage format:	XML, PDF, Mindray
Report Formats:	3x4+1R, 3x4+3R, 6x2, 6x2+1R, 12x1
Additional Functions	
	Preview and review full ECG report on screen
	Reanalyze ECG automatically after changing patient's demographics
	Touch screen (Optional)
	Connect to external printer directly
	Upload XML reports through FTP protocol
	Support barcode reader
	Wifi (Optional)
	USB flash drive storage of PDF and XML outputs
Accessories	
	ECG patient cable with banana plugs, Limb Clamps, Chest Bulbs (IEC/AHA)
	ECG cable with Electrode clips (IEC/AHA)
	ECG patient cable with banana plugs, banana to tab lead adapters, ECG tab electrodes
	Country-specific power cords
	Z-fold thermal paper (Letter or A4)
	Trolley
	Laser printer
	Barcode reader

Mindray Building, Keji 12th Road South,
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com www.mindray.com

Mindray is listed on the NYSE under the symbol "MR"

mindray | healthcare with heart BeneHeart are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD.
©2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.
P/N:ENG-BeneHeart R12-210285x2P-20131112



mindray